

Anti-PRTN3 antibody



Catalog Number: 175103

Product name

Anti-PRTN3 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to PRTN3

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC, IHC-P

Dilutions

WB: 1:500

ICC: 1:50-1:200

IHC-P: 1:100-1:500

Validations




Fig1: Western blot analysis of PRTN3 on U937 using anti-PRTN3 antibody at 1/500 dilution.





Fig2: Western blot analysis of PRTN3 on mouse marrow lysate using anti-PRTN3 antibody at 1/100



dilution.

Fig3: ICC staining PRTN3 in 293T cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig4: ICC staining PRTN3 in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig5: ICC staining PRTN3 in SH-SY5Y cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde,

permeabilised with 0.25% Triton X100/PBS.





Fig6: Immunohistochemical analysis of paraffin-embedded human lung cancer tissue using anti-PRTN3 antibody. Counter stained with



hematoxylin.

Fig7: Immunohistochemical analysis of paraffin-embedded human spleen tissue using anti-PRTN3 antibody. Counter stained with hematoxylin.