

Anti-fezf2 antibody



Catalog Number: 175116

Product name

Anti-fezf2 antibody

Specificity

Zebrafish

Antibody description

Rabbit polyclonal antibody to fezf2

Preparation

This antigen of this antibody was recombinant protein within human fezf2 aa 211-385.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:1000

Validations




Fig1: Western blot analysis of FEZF2 on Zebrafish tissue lysates using anti-FEZF2 antibody at 1/500 dilution.