

Anti-CTDSPL antibody



Catalog Number: 175133

Product name

Anti-CTDSPL antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to CTDSPL

Preparation

This antigen of this antibody was recombinant protein within human ctdspl aa 50-300.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol. Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, FC

Dilutions

WB: 1:500

FC: 1:50-1:100

Validations





Fig1: Western blot analysis of CTDSPL on HepG2 cell lysate using anti-CTDSPL antibody at 1/500



dilution.

Fig2: Flow cytometric analysis of HepG2 cells with CFTR antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti-rabbit IgG was used as the secondary antibody.