

Anti-MICU1 antibody



Catalog Number: 175139

Product name

Anti-MICU1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to MICU1

Preparation

This antigen of this antibody was recombinant protein within human cbara1 aa 14-154.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-1:2,000

Validations




Fig1: Western blot analysis of CBARA1 on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody was used at a 1:500 dilution in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature.; Positive control.; Lane 1: Mouse kidney tissue lysate; Lane 2: Rat brain tissue lysate; Lane 3: HeLa cell lysate