

Anti-AFDN antibody



Catalog Number: 175195

Product name

Anti-AFDN antibody

Specificity

Human, Mouse, Rat, Dog, Pig, Cow, Horse

Antibody description

Rabbit polyclonal antibody to AFDN

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human af6 341-440/1816

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB:1:500-2000

IHC-P:1:400-800

Validations





Fig1: Sample:; Hippocampus (Mouse) Lysate at 40 ug; Primary: Anti-Afadin at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 136 kD;



Observed band size: 135/138 kD

Fig2: Sample:; Skin (Rat) Lysate at 40 ug; Primary: Anti-Afadin at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 136 kD; Observed band size: 135/138 kD




Fig3: Paraformaldehyde-fixed, paraffin embedded (mouse brain tissue); Antigen retrieval by boiling

Anti-AFDN antibody



Catalog Number: 175195

in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation

with (Afadin) Polyclonal Antibody, Unconjugated at 1:400 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.