

General Information

Product name:

Anti-PAEP antibody

Immunogen:

Applications:

WB

Specificity:

Cow, Goat

Antibody description:

Rabbit polyclonal antibody to PAEP

Properties

Preparation:

This antigen of this antibody was beta-lactoglobulin

Formulation:

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage:

Store at -20°C for one year. Avoid repeated

freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality:

Polyclonal

Ig Type:IgG

Applications

WB:1:500-2000

Validations




Fig1: Sample: Goat milk at 30ug;; Primary: Anti-Beta-lactoglobulin at 1:300 dilution;; Secondary: HRP conjugated Goat Anti-Rabbit IgG(bs-0295G-HRP) at 1: 3000 dilution;; Predicted band size : 18kD; Observed band size : 18kD