

Anti-C6orf15 antibody



Catalog Number: 175275

Product name

Anti-C6orf15 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to C6orf15

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human c6orf15 21-120/325

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample: MDA-MB-231 Cell (Human) Lysate at 30 ug; Primary: Anti-C6orf15 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 32kD; Observed band size: 34kD