

Anti-CHST3 antibody

Catalog Number: 175329

Product name

Anti-CHST3 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to CHST3

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human chst3 21-120/479

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB:1:500-2000

IHC-P:1:400-800

Validations





Fig1: Sample: A549 (human) Cell Lysate at 40 ug; Primary: Anti-CHST at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 55 kD; Observed



band size: 55 kD

Fig2: Sample: Spleen (Mouse) Lysate at 40 ug; Primary: Anti-CHST3 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 55 kD; Observed band size: 60 kD




Fig3: Paraformaldehyde-fixed, paraffin embedded (mouse testis tissue); Antigen retrieval by boiling

Anti-CHST3 antibody



Catalog Number: 175329

in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation

with (CHST3) Polyclonal Antibody, Unconjugated at 1:400 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.