

# Anti-CYP4Z2P antibody



Catalog Number: 175362

## Product name

Anti-CYP4Z2P antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to CYP4Z2P

## Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human cyp4z2p 51-150/340

## Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

## Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB:1:500-2000

## Validations




Fig1: Sample:293T (Human)Cell Lysate at 40 ug; Primary: Anti-CYP4Z2P at 1/300 dilution; Secondary: IRDye800CW Goat Anti-RabbitIgG at 1/20000 dilution; Predicted band size: 40kD; Observed band size: 40kD