

Anti-DLGAP4 antibody



Catalog Number: 175378

Product name

Anti-DLGAP4 antibody

Specificity

Human, Mouse, Rat, Dog, Pig, Cow, Horse, Rabbit, Sheep, Guinea Pig

Antibody description

Rabbit polyclonal antibody to DLGAP4

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human sapap4/dlgap4 261-370/992

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; Muscle (Mouse) Lysate at 40 ug; Primary: Anti- DLGAP4 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 108 kD; Observed band size: 108 kD