

Anti-NRROS antibody



Catalog Number: 175560

Product name

Anti-NRROS antibody

Specificity

Human, Mouse, Rat, Rabbit

Antibody description

Rabbit polyclonal antibody to NRROS

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human *lrrc33* 501-600/692

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF

Dilutions

WB:1:500-2000

IF:1:100-500

ICC/IF:1:100-500

Validations




Fig1: Sample:; Bone(Rat) Cell Lysate at 40 ug; Thymus(Rat) Cell Lysate at 40 ug; Primary: Anti-LRRC33 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 74 kD; Observed band size: 74 kD