

Anti-Myod1 antibody



Catalog Number: 175586

Product name

Anti-Myod1 antibody

Specificity

Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Rabbit, Sheep

Antibody description

Rabbit polyclonal antibody to Myod1

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from mouse myod1 51-150/318

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; Muscle (Rat) Lysate at 40 ug; Muscle (Mouse) Lysate at 40 ug; Primary: Anti-MyoD1/Myf3 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 35 kD; Observed band size: 36 kD