

Anti-OVAL antibody

Catalog Number: 175625

Product name

Anti-OVAL antibody

Specificity

Chicken

Antibody description

Rabbit polyclonal antibody to OVAL

Preparation

This antigen of this antibody was ovalbumin [hen egg white]

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; Ovalbumin protein at 5 ng; Ovalbumin protein at 5 ng; Primary: Anti-Ovalbumin at 1/500 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 43 kD; Observed band size: 50 kD