

Anti-TRIM34 antibody



Catalog Number: 175668

Product name

Anti-TRIM34 antibody

Specificity

Human, Mouse, Rat, Pig

Antibody description

Rabbit polyclonal antibody to TRIM34

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human rnf21 301-400/488

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; Lung (Mouse) Lysate at 40 ug; Placenta (Mouse) Lysate at 40 ug; Primary: Anti-RNF21 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 57kD; Observed band size:




Fig2: Sample:; Lung (Mouse) Lysate at 40 ug; Primary: Anti- RNF21 at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 57kD; Observed band size: 55 kD