

Anti-LOC102724334 antibody



Catalog Number: 175803

Product name

Anti-LOC102724334 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to LOC102724334

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human histone h2b around the acetylated site of lys15. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of 3T3 cells using Acetyl-Histone H2B (K15) Polyclonal Antibody diluted at 1: 1000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000