

# Anti-H1-1 antibody



Catalog Number: 175805

## Product name

Polyclonal

Anti-H1-1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to H1-1

## Dilutions

## Preparation

WB: 1:500-1:1000

This antigen of this antibody was synthetic peptide of histone h1 (di methyl lys25)

## Validations

## Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1.; Western blot analysis of 1) HeLa, 2) 3T3, diluted at 1:2000.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000