

# Anti-ZBTB45 antibody

Catalog Number: 175807

## Product name

Anti-ZBTB45 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to ZBTB45

## Preparation

This antigen of this antibody was recombinant protein of zbtb45

## Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:1000

IHC-P: 1:200-1:500

## Validations




Fig1.; Western blot analysis of 1) Helix, 2) Jurkat, 3) Mouse Brain, 4) Mouse Kidney, 5) Rat Brain using ZBTB45 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at




Fig2.; Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using ZBTB45 Polyclonal Antibody.