

# Anti-BEST2 antibody

Catalog Number: 175810

## Product name

Anti-BEST2 antibody

## Specificity

Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to BEST2

## Preparation

This antigen of this antibody was synthetic peptide of best2

## Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:1000

IHC-P: 1:200-1:500

## Validations




Fig1:: Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue using BEST2 Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using BEST2 Polyclonal Antibody.