

# Anti-IGFN1 antibody

Catalog Number: 175812

## Product name

Anti-IGFN1 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to IGFN1

## Preparation

This antigen of this antibody was recombinant protein of igfn1

## Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:2000-1:5000

IHC-P: 1:200-1:500

## Validations




Fig1:: Western blot analysis of 1) Pig Skeletal MuscleTissue, 2) Rat Skeletal Muscle Tissue, 3) Rat Skeletal Muscle Tissue using IGFN1 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded Mouse Skeletal MuscleTissue using IGFN1 Polyclonal Antibody.