

Anti-ADGRB1 antibody

Catalog Number: 175813

Product name

Anti-ADGRB1 antibody

Specificity

Human, Rat, Mouse

Antibody description

Rabbit polyclonal antibody to ADGRB1

Preparation

This antigen of this antibody was synthetic peptide of bail

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1::; Western blot analysis of 1) 293T, 2) Mouse Brain Tissue, 3) Rat Brain Tissue with BAI1 Rabbit pAb diluted at 1:2,000.




Fig2::; Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using BAI1 Rabbit pAb diluted at 1:200.