

Anti-CACNA1I antibody



Catalog Number: 175817

Product name

Anti-CACNA1I antibody

Specificity

Human, Rat, Mouse

Antibody description

Rabbit polyclonal antibody to CACNA1I

Preparation

This antigen of this antibody was synthetic peptide of cav3.3

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

IHC-P

Dilutions

IHC-P: 1:100-200

Validations



Fig1.; Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Cav3.3Rabbit pAb diluted at 1:200.