

Anti-KCNN2 antibody

Catalog Number: 175819

Product name

Anti-KCNN2 antibody

Specificity

Human, Rat, Mouse

Antibody description

Rabbit polyclonal antibody to KCNN2

Preparation

This antigen of this antibody was synthetic peptide of kcnn2 (sk2)

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:1000-1:2000

IHC-P: 1:100-200

Validations




Fig1:: Western blot analysis of 1) Rat BrainTissue, 2) Mouse Brain Tissue, 3) HepG2 with KCNN2(SK2) Rabbit pAb diluted at 1:2,000.




Fig2:: Immunohistochemical analysis of paraffin-embedded Human BrainTissue using KCNN2(SK2) Rabbit pAb diluted at 1:200.