

# Anti-KCNG3 antibody

Catalog Number: 175820

---

## Product name

Anti-KCNG3 antibody

## Specificity

Human, Rat, Mouse

## Antibody description

Rabbit polyclonal antibody to KCNG3

## Preparation

This antigen of this antibody was synthetic peptide of kv10.1

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

IHC-P

## Dilutions

IHC-P: 1:100-200

## Validations




Fig1.; Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Kv10.1 Rabbit pAb diluted at 1:200.