

# Anti-DBX1 antibody



Catalog Number: 175823

## Product name

Anti-DBX1 antibody

## Specificity

Human, Rat, Mouse

## Antibody description

Rabbit polyclonal antibody to DBX1

## Preparation

This antigen of this antibody was synthesized peptide derived from part region of human protein

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1:: Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°over night.