

Anti-DBX1 antibody



Catalog Number: 175823

Product name

Anti-DBX1 antibody

Specificity

Human, Rat, Mouse

Antibody description

Rabbit polyclonal antibody to DBX1

Preparation

This antigen of this antibody was synthesized peptide derived from part region of human protein

Formulation

Liquid, Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°over night.