

Anti-AAK1 antibody



Catalog Number: 175830

Product name

Anti-AAK1 antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Rabbit polyclonal antibody to AAK1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human aak1. aa range:271-320

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations



Fig1.; Western Blot analysis of various cells using AAK1 Polyclonal Antibody