

# Anti-ACTG2 antibody



Catalog Number: 175840

## Product name

Anti-ACTG2 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to ACTG2

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human actin-gamma2. aa range:1-50

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IF

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

## Validations




Fig1.;; Western Blot analysis of various cells using Actin  $\alpha$ 3 Polyclonal Antibody diluted at 1: 2000