

Anti-ARSK antibody

Catalog Number: 175864

Product name

Anti-ARSK antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to ARSK

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human arsk. aa range:481-530

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of K562 cells using Arylsulfatase K Polyclonal Antibody