

Anti-CDH7 antibody



Catalog Number: 175880

Product name

Anti-CDH7 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to CDH7

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cdh7. aa range:651-700

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations



Fig1.; Western Blot analysis of various cells using Cadherin-7 Polyclonal Antibody diluted at 1: 1000