

Anti-CAPN11 antibody

Catalog Number: 175884

Product name

Anti-CAPN11 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CAPN11

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human capn11. aa range:401-450

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations

Fig1.; Western Blot analysis of various cells using Calpain 11 Polyclonal Antibody