

# Anti-CNGB1 antibody



Catalog Number: 175922

## Product name

Polyclonal

Anti-CNGB1 antibody

## Ig Type

## Specificity

IgG

Human, Mouse, Rat

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to CNGB1

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cngb1. aa range:571-620

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

## Validations




Fig1:; Western Blot analysis of various cells using CNG-1β Polyclonal Antibody