

Anti-CYB561D2 antibody



Catalog Number: 175933

Product name

Anti-CYB561D2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CYB561D2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human c56d2. aa range:55-104

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of Jurkat cells using CYB561D2 Polyclonal Antibody