

Anti-CNNM2 antibody

Catalog Number: 175935

Product name

Anti-CNNM2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to CNNM2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human *cnnm2*. aa range:571-620

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000


Validations

Fig1.; Western Blot analysis of various cells using Cyclin M2 Polyclonal Antibody