

# Anti-CYP2A7 antibody



Catalog Number: 175938

## Product name

Polyclonal

Anti-CYP2A7 antibody

## Ig Type

## Specificity

Rabbit IgG

Human

## Applications

## Antibody description

WB, IHC-P

Rabbit polyclonal antibody to CYP2A7

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was synthesized peptide derived from the c-terminal region of human cyp2a7.

IHC-P: 1:100-1:300

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1.; Western Blot analysis of various cells using CYP2A7 Polyclonal Antibody diluted at 1: 2000