

# Anti-CYP4Z1 antibody

Catalog Number: 175942

## Product name

Anti-CYP4Z1 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to CYP4Z1

## Preparation

This antigen of this antibody was synthesized peptide derived from the internal region of human cyp4z1/2.

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

IHC-P: 1:100-1:300

## Validations




Fig1.; Western Blot analysis of various cells using CYP4Z1/2 Polyclonal Antibody