

Anti-DCLK2 antibody

Catalog Number: 175945

Product name

Anti-DCLK2 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to DCLK2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human dclk2. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations



Fig1.; Western Blot analysis of various cells using DCAMKL2 Polyclonal Antibody