

Anti-POLA1 antibody

Catalog Number: 175952

Product name

Polyclonal

Anti-POLA1 antibody

Ig Type**Specificity**

Rabbit IgG

Human

Applications**Antibody description**

WB

Rabbit polyclonal antibody to POLA1

Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human dna polymerase alpha. aa range:81-130

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Clonality

Fig1.; Western Blot analysis of various cells using DNA pol α Polyclonal Antibody diluted at 1: 1000