

Anti-MEGF8 antibody



Catalog Number: 175956

Product name

Anti-MEGF8 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to MEGF8

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human megf8. aa range:1243-1292

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using EGFL4 Polyclonal Antibody