

Anti-EOMES antibody



Catalog Number: 175963

Product name

Anti-EOMES antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to EOMES

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human eomes. aa range:59-108

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using EOMES Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).