

Anti-EPS8L1 antibody



Catalog Number: 175967

Product name

Anti-EPS8L1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to EPS8L1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human es8l1. aa range:411-460

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using Eps8L1 Polyclonal Antibody