

Anti-FOXD4L1 antibody



Catalog Number: 175983

Product name

Anti-FOXD4L1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to FOXD4L1

Preparation

This antigen of this antibody was synthesized peptide derived from foxd4l1 . at aa range: 320-400

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of HuvEc cells using FoxD4L1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).