

Anti-FOXN4 antibody



Catalog Number: 175986

Product name

Anti-FOXN4 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to FOXN4

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human foxn4. aa range:411-460

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of LOVO cells using FoxN4 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotec, MN, USA).