

Anti-GPR101 antibody

Catalog Number: 175992

Product name

Anti-GPR101 antibody

Specificity

Human, Mouse, Monkey

Antibody description

Rabbit polyclonal antibody to GPR101

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gpr101. aa range:451-500

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using GPR101 Polyclonal Antibody diluted at 1: 500