

Anti-SLC52A1 antibody

Catalog Number: 176003

Product name

Anti-SLC52A1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to SLC52A1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human pcvr2. aa range:235-284

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of 293T cells using GPR172B Polyclonal Antibody diluted at 1: 1000