

Anti-FFAR1 antibody



Catalog Number: 176010

Product name

Polyclonal

Anti-FFAR1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Monkey

Applications

Antibody description

WB, IF

Rabbit polyclonal antibody to FFAR1

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ffar1. aa range:185-234

IF: 1:100-1:500

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Validations

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1:; Western Blot analysis of COS7 cells using GPR40 Polyclonal Antibody diluted at 1: 500