

Anti-FFAR2 antibody



Catalog Number: 176011

General Information

Product name:

Anti-FFAR2 antibody

Immunogen:

Applications:

IF

Specificity:

Human

Antibody description:

Rabbit polyclonal antibody to FFAR2

Properties

Preparation:

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ffar2. aa range:271-320

Formulation:

Liquid, Liquid in PBS containing 50% glycerol,

0.5% BSA and 0.02% sodium azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality:

Polyclonal

Ig Type: IgG

Applications

IF: 1:100-1:500

Validations




Fig1.; Immunofluorescence analysis of COS7 cells, using FFAR2 Antibody. The picture on the right is blocked with the synthesized peptide.