

Anti-SUCNR1 antibody



Catalog Number: 176013

Product name

Anti-SUCNR1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to SUCNR1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human sucnr1. aa range:100-149

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western blot analysis of lysates from HUVEC and MCF-7 cells, using SUCNR1 Antibody. The lane on the right is blocked with the synthesized peptide.