

Anti-GRK7 antibody

Catalog Number: 176017

Product name

Anti-GRK7 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to GRK7

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human grk7. aa range:371-420

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of various cells using GRK 7 Polyclonal Antibody