

Anti-MAPK8IP3 antibody

Catalog Number: 176037

Product name

Anti-MAPK8IP3 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to MAPK8IP3

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human jip3. aa range:621-670

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations (kD)




Fig1.; Western Blot analysis of various cells using JIP-3 Polyclonal Antibody diluted at 1: 1000