

# Anti-BRI3BP antibody



Catalog Number: 176045

## Product name

Anti-BRI3BP antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to BRI3BP

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human bri3b. aa range:1-50

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IF

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

## Validations




Fig1:: Western blot analysis of lysates from RAW264.7 cells, using BRI3B Antibody. The lane on the right is blocked with the synthesized peptide.




Fig2:: Immunohistochemistry analysis of paraffin-embedded human brain tissue, using BRI3B Antibody. The picture on the right is blocked with the synthesized peptide.